#### REPUBLIK INDONESIA KEMENTERIAN HUKUM DAN HAK ASASI MANUSIA

# SURAT PENCATATAN CIPTAAN

Dalam rangka pelindungan ciptaan di bidang ilmu pengetahuan, seni dan sastra berdasarkan Undang-Undang Nomor 28 Tahun 2014 tentang Hak Cipta, dengan ini menerangkan:

Indonesia

Indonesia

:

Buku Panduan/Petunjuk

29 April 2024, di Jakarta Barat

Januari tahun berikutnya.

: EC00202433771, 29 April 2024

Dr. dr. Siufui Hendrawan, M.Biomed

Dr. dr. Siufui Hendrawan, M.Biomed

Jl. Letjen S. Parman No. 1, Universitas Tarumanagara, Fakultas Kedokteran, Grogol Petamburan, Jakarta Barat, DKI Jakarta, 11440

Jl. Letjen S. Parman No. 1, Universitas Tarumanagara, Fakultas

Kedokteran, Grogol Petamburan, Jakarta Barat, Dki Jakarta 11440

**Preparation Of Collagen-Coated PLLA Matrix Seeded With** 

Berlaku selama hidup Pencipta dan terus berlangsung selama 70 (tujuh puluh) tahun setelah Pencipta meninggal dunia, terhitung mulai tanggal 1

Fibroblast And Supplemented With Secretome

Nomor dan tanggal permohonan

#### Pencipta

Nama

Alamat

Kewarganegaraan

Pemegang Hak Cipta

Nama

Alamat

Kewarganegaraan

Jenis Ciptaan

Judul Ciptaan

Tanggal dan tempat diumumkan untuk pertama kali di wilayah Indonesia atau di luar wilayah Indonesia

Jangka waktu pelindungan

Nomor pencatatan

adalah benar berdasarkan keterangan yang diberikan oleh Pemohon. Surat Pencatatan Hak Cipta atau produk Hak terkait ini sesuai dengan Pasal 72 Undang-Undang Nomor 28 Tahun 2014 tentang Hak Cipta.

000609127



a.n. MENTERI HUKUM DAN HAK ASASI MANUSIA DIREKTUR JENDERAL KEKAYAAN INTELEKTUAL u.b T<sub>1</sub>I

 $\mathbf{V}$ 

Direktur Hak Cipta dan Desain Industri

IGNATIUS MT SILALAHI NIP. 196812301996031001

Disclaimer:

Dalam hal pemohon memberikan keterangan tidak sesuai dengan surat pernyataan, Menteri berwenang untuk mencabut surat pencatatan permohonan.

V

0	TARUMANAGARA HUMAN CELL TECHNOLOGY L A B O R A T O R Y THCT Lab	TARUMANAGARA HUMAN CELL TECHNOLOGY (THCT) LABORATORY FACULTY OF MEDICINE, TARUMANAGARA UNIVERSITY Jl. Letjen S. Parman No. 1; Jakarta 11440 INDONESIA Phone. +62 21 5696 3254, Fax. +62 21 5696 7325, Email. thctlab11@gmail.com	
		STANDARD OPERATING PROCEDURE (SOP)	Document Num.: S-031 Release date: 25 April 2024
			Effective date: 25 April 2024
		Preparation of Collagen-Coated PLLA Matrix Seeded with Fibroblast and Supplemented with Secretome	Version: 1
			First draft: THCT
			Approved by:

# 1. Preparation of Culture Medium Supplemented with Secretome (DMEM+Sec)

- Dulbecco's Modified Eagle's Medium (DMEM; Sigma) complete medium is prepared by adding 10% Fetal Bovine Serum (FBS; Gibco), 1% sodium pyruvate (Sigma), 1% 4-(2-Hydroxyethyl) piperazine-1-ethanesulfonic acid (HEPES; Sigma), and 1% Antibioticantimycotic solution (Ab/Am; Sigma).
- For secretome supplementation (DMEM+Sec), DMEM complete medium is added with 1% secretome derived from human umbilical cord mesenchymal stem cells (hUC-MSCs).

Note: hUC-MSCs secretome is obtained from the protocol of HKI No. 000284674.

## 2. Production of Ready to Use Matrix for Hernia Repair

Collagen-coated poly-L-lactide (PLLA; Resomer L206S, Evonik) matrix is fabricated through salt-leaching method, following the procedure from Hendrawan *et al.* (2020)<sup>1</sup> (Fig 1).

<u>Note:</u> The matrix is cut into a rectangular shape of 2 x 1 cm with blunt edges.

- $\circ$  One day before the hernia repair and implant surgery, each sterile matrix is seeded with 400 μL of fibroblast suspension, containing 2x10<sup>6</sup> cells.
- After cell seeding, 800 uL of DMEM+Sec is added into the surrounding of matrix and further incubated for 1 hour at 37°C in a 5% CO<sub>2</sub> incubator.
- After 1 hour of incubation, another 1.2 mL of DMEM+Sec medium is added into each well and then incubated for 24 hours.
- $\circ$  Prior to hernia repair surgery, the matrix is transferred into a new 6 well-plate containing 800  $\mu$ L of DMEM complete medium. The matrix is ready for use.



Figure 1 Rectangular collagen-coated PLLA matrix with 2 cm length (A) and 1 cm width (B).

## 3. Hernia Repair and Implant Treatment Procedure

Method is modified from Hendrawan et al. 2024<sup>2</sup>.

- Male Wistar rat weighing 180-200 gr is anesthetized using 10% Ketamine (40-80 mg/kg BW) and 2% Xylazine (5-10 mg/kg BW).
- The hair on ventral abdominal side is removed using hair clipper and blade razor. The skin surface is disinfected with 70% ethanol and povidone-iodine.
- A midline laparotomy of ±2 cm in length is done along the *linea alba* to create hernia in the animal model (Fig 2A).
  Note: The abdominal muscle incision must be done gently without rupturing the

<u>Note:</u> The abdominal muscle incision must be done gently without rupturing the peritoneum membrane.

- At the site of laparotomy, either side of the muscle is carefully separated from the peritoneum membrane to give a space for the matrix (Fig 2B).
- The implant is gently placed in sublay (in between of abdominal muscle and peritoneum membrane), without breaking the implant (Fig 2C).
- Abdominal muscle and skin are closed with a continuous suture using 4-0 absorbable vicryl (Ethicon W9386) (Fig 2D).
- Skin surgical site is applied with Gentamicin ointment and covered with Hypafix.
- Rats are kept until 2 months, for further tensile strength testing and histopathological staining of the abdominal muscle.



**Figure 2** Procedure for hernia repair and implant treatment in rat model: midline laparotomy of 1.8 cm is created on the abdominal muscle (A), the peritoneum membrane is gently separated from the muscle (B), the implant is inserted sublay between the muscle and peritoneum (C), abdominal muscle and skin are sutured (D).

# 4. Euthanasia and Tissue Collection

- The rat is anesthetized using 10% Ketamine (80 mg/kg BW) and 2% Xylazine (10 mg/kg BW).
- Under deep anesthesia, rat is euthanized through intracardiac injection of sodium pentobarbital (200 mg/kg BW).
- Ventral abdominal side of the rat is shaved using a hair clipper.
- The ventral abdominal skin is separated from the muscle part and discarded (Fig 3A).
- The abdominal muscle on the area of midline laparotomy and implant is obtained according to the corresponding test:
  - <u>For tensile strength analysis:</u> a 6x4 cm of muscle is collected and submerged in a physiological salt solution for immediate tensile pull testing (Fig 3B).
  - For histopathological staining: a 3x2 cm of muscle is collected and submerged in a 4% paraformaldehyde solution to be made into a paraffin block (Fig 3C).



**Figure 3** Sample collection and preparation of abdominal muscle: separation of abdominal skin from the muscle (A), collection of abdominal muscle for tensile strength analysis (B) and histopathological staining (C).

# 5. Tensile Strength Analysis of Abdominal Muscle at The Implanted Area

Method is modified from Melman et al. (2011)<sup>3</sup> and Holmdahl et al. (2019)<sup>4</sup>.

- Tensile strength is measured using a universal testing machine with a dynamic force processor (Unipulse F381A).
- Obtained tissue is stretched and fixed into the upper and lower grips, with the direction of midline laparotomy is positioned laterally (Fig 4).
  <u>Note:</u> The lower grip should hold the implant, while the upper grip holds the tissue near the midline laparotomy. Make sure that the tissue is all stretched without any parts folded.
- A uniform speed of 400 rpm is applied and the pulling test is performed until the tissue ruptures.
- $\circ$  The maximum force borne by each specimen is recorded (in Newton).



Figure 4 Tensile strength testing setup.

## 6. Collagen I:III Ratio Confirmation

To ascertain the effect of treatment for a stronger muscle development, collagen observation using Sirius Red solution is a widely-performed histology staining. Particularly through this staining, deposition of collagen type I and type III could be differentially observed using a polarized microscope. Collagen type I is formed as a thicker and stronger bundle, while collagen type III deposition is associated with a weaker muscle development. Hence, the ratio of collagen I:III could depict the strength of hernia repair<sup>5</sup>.

- To distinguish the deposition of collagen type I from the type III, the stained specimen is observed under a polarized light using a light microscope (Olympus CX23) fixed with a custom-made analyzer and polarizer.
- Pictures are taken in 10x field views for further determination of collagen I:III ratio.
- The determination of ratio collagen I:III is quantified using ImageJ software (National Institute of Health)<sup>6</sup>.

## References

- Hendrawan, S. *et al.* Evaluation of 3D PLLA scaffolds coated with nano-thick collagen as carrier for hepatocytes. *J. Biomed. Mater. Res. - Part B Appl. Biomater.* **109**, 723–732 (2021).
- 2. Hendrawan, S. *et al.* Fibroblast matrix implants—a better alternative for incisional hernia repair? *Biomed. Mater.* **19**, 035033 (2024).
- 3. Melman, L. *et al.* Early biocompatibility of crosslinked and non-crosslinked biologic meshes in a porcine model of ventral hernia repair. *Hernia* **15**, 157–164 (2011).
- Holmdahl, V., Backman, O., Gunnarsson, U. & Strigård, K. The Tensile Strength of Full-Thickness Skin : A Laboratory Study Prior to Its Use as Reinforcement in Parastomal Hernia Repair. 6, 2–7 (2019).
- 5. Benito-Martínez, S. *et al.* Self-adhesive hydrogel meshes reduce tissue incorporation and mechanical behavior versus microgrips self-fixation: a preclinical study. *Hernia* **26**, 543–555 (2022).
- 6. Fuentes, R., Saravia, D., Arias, A., Prieto, R. & Dias, F. Mandibular dental implant placement using demineralized bone matrix (DBM). *Biomed. Res.* 28, 2656–2660 (2017).