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 <p>TARUMANAGARA HUMAN CELL TECHNOLOGY LABORATORY THCT Lab</p>	TARUMANAGARA HUMAN CELL TECHNOLOGY (THCT) LABORATORY FACULTY OF MEDICINE, TARUMANAGARA UNIVERSITY Jl. Letjen S. Parman No. 1; Jakarta 11440 INDONESIA Phone. +62 21 5696 3254, Fax. +62 21 5696 7325, Email. thctlab11@gmail.com	
	STANDARD OPERATING PROCEDURE	Document Num.: S-032
	(SOP)	Release date: 25 Sep 2024
		Effective date: 25 Sep 2024
	Protein Isolation from Rat Liver Tissue (Cryopreserved Organ)	Version: 1
		First draft: THCT
		Approved by:

1. Liver Tissue Procurement and Storage

- During necropsy, the liver organ is cut into smaller tissue sections.
- Part of liver section is weighed and separated into a new sterile tube.

Note: For protein isolation, usually only ± 0.1 g of liver section is required for each 0.5–1 mL of lysis buffer cocktail.

- The organ is immediately placed in a cool box containing ice prior to storage.
- For long-term storage, the organ could be stored in deep freezer at -80°C.

2. Tissue Pulverization^{1,2}

All process should be performed on ice block or at 4°C

- RIPA/PI cocktail is prepared by adding 10 µL protease inhibitor (PI; MCE Protease Inhibitor Cocktail #HY-K0010) to every 1 mL of RIPA lysis buffer (Elabscience #E-BC-R327) in a sterile microcentrifuge tube.

Note: Pre-calculate the volume of buffer cocktail required for all samples.

Each 0.1 g of tissue sample could be added with 0.5 – 1 mL lysis buffer (1:5 – 1:10 ratio), adjusted according to preferred experimental requirements. The RIPA/PI cocktail should be prepared freshly, right before adding into the sample. Keep the cocktail on cooling block or place it in -20°C while waiting.

- Mortar, pestle, and spatula are precooled with liquid nitrogen (LN) by carefully pouring the LN into the mortar and swirling the LN slowly using the pestle and spatula, until the LN evaporates thoroughly. This step is performed 2-3 times.

Note: Clean all the utensils thoroughly with 70% alcohol beforehand.

- Adequate amount of LN is poured into the cold mortar, followed by soaking the frozen tissue sample in LN.
- The tissue was ground immediately until become a fine powder; LN could be added accordingly if needed.

Note: Grinding the tissue in LN should be done carefully but also as quick as possible to prevent the temperature rising.

- RIPA/PI cocktail of 0.5 – 1 mL is added into the tissue powder, continued by grinding and mixing it again.
- The powder is scraped from all over the mortar using a cold spatula and transferred into a sterile microcentrifuge tube.

Note: Keep the whole process in a cold condition. Place the microcentrifuge tube in a cooling block or inside a cool box containing ice.

3. Protein Isolation

- Protein sample is incubated on ice or at 4°C for 20 minutes
- Sample is centrifuged at 12,000 rpm, 4°C for 10 minutes.
- The supernatant is collected and transferred into a new sterile tube.

Note: Carefully collect the supernatant without disturbing the pellet.

- Sample is then subjected to protein concentration or purity measurement. If the concentration is too high, PBS solution may be added to dilute the sample.

References

1. Elabscience Bionovation Inc. RIPA Lysis Buffer (Strong). *Product Catalog 1–2*.
2. Menggensilimu *et al.* Anti-liver fibrosis effect of total flavonoids from Scabiosa comosa Fisch. ex Roem. et Schult. on liver fibrosis in rat models and its proteomics analysis. *Ann. Palliat. Med.* **9**, 272–285 (2020).